

# Fiber Optic Tweezers and Optical Stretcher: biophotonic tools for single cell studies

F. Bragheri, L. Ferrara , P. Minzioni, I. Cristiani

*Electronics Department, University of Pavia, Italy*

francesca.bragheri@unipv.it



Regione Lombardia

Pavia, 17/11/2009

Francesca Bragheri



# Motivation

In molecular and cellular biology studies performed at a single cell level have revealed extremely **heterogeneous behaviours** and provided unpredictable information

The exploitation of **optical forces**, generated by the interaction between electromagnetic waves and the sample, represents an accurate, non-invasive and gentle manipulation technique for individual cell studies

**Bio-photonics devices** based on optical forces:

- can trap and move single cells **without physical contact**
- can be combined with diagnostic methods
- can be used for many applications including, cell sorting, bacterial adhesion forces and membrane interactions analysis and investigation of the **viscoelastic properties**

**Aim: realization of miniaturized biophotonics tools through the use of optical fibers and integrated optics**

# Outline


- Biophotonic tools based on optical forces:  
Tweezers' working principle
- Biophotonic tools based on optical forces:  
Stretcher's working principle
- Optical Stretcher Capabilities and Experimental results
- Conclusion

# Outline

- Biophotonic tools based on optical forces:  
Tweezers' working principle
- Biophotonic tools based on optical forces:  
Stretcher's working principle
- Optical Stretcher Capabilities and Experimental results
- Conclusion



# Light can move the objects!



Micky mouse is pushing the mirror


If light changes its momentum due to the scattering from an object, at the same time it is **applying a force**




# Single-beam Optical Trap

Optical trap: stable equilibrium point of optical forces acting on a particle

The trap is created by a single laser beam tightly focused: large intensity gradients in both the axial and transverse directions



# Standard vs Fiber Optic Tweezers



Standard Optical Tweezers



# Standard vs Fiber Optic Tweezers

## NEWS & VIEWS

### OPTICAL TWEEZERS


## All-fibre design

Optical tweezers enable precise, controlled and non-contact manipulation of small biological specimens. Rather than using a bulky microscope, it is now possible to create optical tweezers at the end of a fibre probe.

#### Miles Padgett

is in the Department of Physics and Astronomy,  
University of Glasgow, Glasgow, G12 8QQ, UK.  
e-mail: m.padgett@physics.gla.ac.uk

In the 20 years since they were pioneered by Arthur Ashkin<sup>1</sup> and his co-workers, optical tweezers have become a mainstream research tool in both the biological and physical sciences. Most optical tweezers rely on the force acting on a dielectric particle that arises whenever there is a large gradient in the intensity of the optical field. Usually this gradient is achieved by the strong focusing of a conventional laser beam using a high-magnification microscope objective. Indeed, most optical tweezers are themselves based around high-quality microscopes where the same objective lens is used for both focusing the trapping beam and viewing the trapped sample. Now, Carlo Liberale and colleagues<sup>2</sup> have come up with a different approach, one that enables tweezing using a single optical-fibre probe (page 723).




Images of the all-fibre tweezer<sup>2</sup>. The tweezer lying next to a syringe (left) and a scanning electron microscope image of the tweezer's end showing four internal fibre structures with specially shaped core regions (right). The diameter of the fibre probe is 360  $\mu\text{m}$ .



Our proposal: Fiber Optic Tweezers



# Fiber Optic Tweezers - Nanofabrication



With a Focused Ion Beam (FIB), we performed micro machining, digging only the core regions where the total internal reflection takes place.



Regione Lombardia


In collaboration with Università della Magna Graecia - Catanzaro and IIT - Genova

Pavia, 17/11/2009

Francesca Bragheri



# Fiber Optic Tweezers - Trapping



**C. Liberale et al. "Miniaturized all-fibre probe for three-dimensional optical trapping and manipulation", Nature Photonics 1, 723, 2007**




Regione Lombardia

Pavia, 17/11/2009


Francesca Bragheri



# Fiber Optic Tweezers - Trapping



# Fiber Optic Tweezers - Fluorescence collection




# Outline

- Biophotonic tools based on optical forces:  
Tweezers' working principle
- Biophotonic tools based on optical forces:  
Stretcher's working principle
- Optical Stretcher Capabilities and Experimental results
- Conclusion

# Dual-beam Optical Trap

Optical trap: stable equilibrium point of optical forces acting on a particle

The easiest way to realize a trap:  
exploitation of two counterpropagating laser beams



**A. Ashkin "Acceleration and trapping of particles by radiation pressure" Phys. Rev. Lett. 24, 156, 1970**



Regione Lombardia

Pavia, 17/11/2009


Francesca Bragheri



# Optical Forces

Refraction of light: change in the carried momentum  
Equal and opposite momentum is imparted to the sphere

$$\text{Force: } \frac{dp}{dt}$$




By illuminating a spherical sample through a gaussian beam two forces are produced:

- scattering force that tends to push the particle
- gradient force that tends to pull the particle toward the beam axis



# Optical Trapping: working principle

When two counterpropagating beams are present, the scattering forces are balanced whereas the gradient components are acting in the same direction






# Outline


- Biophotonic tools based on optical forces:  
Tweezers' working principle
- Biophotonic tools based on optical forces:  
Stretcher's working principle
- **Optical Stretcher Capabilities** and Experimental results
- Conclusion

# Trapping


We efficiently trapped different dielectric and biological samples  
(stem cells, yeast cells, red blood cells, white blood cells...)



- polystyrene bead  
( $d \sim 10 \mu\text{m}$ )




- stem cell  
( $d \sim 10 \mu\text{m}$ )



- red blood cell  
( $d \sim 6-10 \mu\text{m}$ )

Collaborations with groups of the internal network of the project

# Trapping




Regione Lombardia

Pavia, 17/11/2009

Francesca Bragheri



# Moving



Regione Lombardia



Pavia, 17/11/2009

Francesca Bragheri



# Dual-beam Optical Trap: Optical Stretcher

An optical stretcher is a device that can induce mechanical deformations on a sample without physical contact



Red Blood Cells

J. Guck et al. "The optical stretcher: a novel laser tool to micromanipulate cells" *Biophys. J.* 81,767, 2001




Regione Lombardia

Pavia, 17/11/2009

Francesca Bragheri



# Stretching - 1



Hypotonic buffer: Swollen Red Blood Cell



# Stretching - 2



Isotonic buffer: Red Blood Cell




# Outline

- Biophotonic tools based on optical forces:  
Tweezers' working principle
- Biophotonic tools based on optical forces:  
Stretcher's working principle
- **Optical Stretcher** Capabilities and **Experimental results**
- Conclusion




# All-fiber experimental setup




Microfluidic circuit  
used to control the sample flow  
between the two fibers



# Integrated experimental setup



# Experimental results - Diabetes



Diabetic patients with complications:  
7 samples  
(50 RBCs each)


Diabetic patients:  
4 samples  
(50 RBCs each)

Healthy patients:  
4 samples  
(50 RBCs each)

Slope of interpolating curves

	Increasing Power	Decreasing Power
Healthy	0.00115861	0.0013566
Diabetics	0.00096231	0.00077057
Diabetics with comp.	0.00091648	0.00088093

# Experimental results – Rare Genetic Disease



In collaboration with Dipartimento di Patologia Umana ed Ereditaria - Prof C. Danesino and coworkers

# Outline

- Biophotonics tools based on optical forces:  
Tweezers' working principles
- Biophotonics tools based on optical forces:  
Stretcher's working principles
- Optical Stretcher Capabilities and Experimental results
- **Conclusion**

# Conclusions


- Optical forces are a powerful tool for **non-invasive cell manipulation**
  - A **Fiber Optic Tweezers** capable of trapping and manipulating single cells has been realized
  - An **Optical Stretcher** apparatus capable of trapping, manipulating and stretching particles has been realized
- Very preliminar results on different diseases shows a slightly **reduced deformability of patients' red blood cells** with respect to healthy controls
- Optical stretching can be combined to optical diagnostic techniques for **cell sorting and analysis**
- We believe that **integrated optics approach** could provide reliable and handy tools for real applications

## Thanks for your attention!


francesca.bragheri@unipv.it

# Optical Stretching: working principle





Single ray impinging on a particle



Two counterpropagating beams




# Deformation analysis





# REALIZZAZIONE DELL'APPARATO


Tecnica per il controllo del flusso



Alzando, con un traslatore micrometrico, di  $\Delta h$  uno dei due tubicini si produce il flusso di una quantità di soluzione pari a  $S\Delta h$ , dove  $S$  è la sezione del tubicino

# Application in biology

Investigation of the viscoelastic properties through the application of intense optical forces, able to cause a significant **deformation of the cytoskeleton** whose degree of deformability is altered by many diseases and provides a **marker of the cell status**



Deformation study on human breast epithelial cells:

- MCF-10 non tumorigenic epithelial cells
- MCF-7 human breast cancer cell non-methastatic
- Mod-MCF-7 methastatic cell







FIGURE 6 Typical examples of the stretching of breast epithelial cells. The images in the left column are taken at an incident light power of 100 mW in each beam, which is sufficient for the trapping of the cells. At an incident light power of 600 mW (right column), the cancerous MCF-7 cells (C and D) deform more than the nonmalignant MCF-10 cells (A and B). The metastatic modMCF-7 cells (E and F) deform the most. The scale bar is 10  $\mu$ m.

- J. Guck et al "Optical Deformability as Inherent Cell Marker for Malignant Transformation and Metastatic Competence" *Biophys. J.* 88, 3689, 2005.







Regione Lombardia

Pavia, 17/11/2009

Francesca Bragheri

